

## Product information

### In.QualiCheck Diabetes GAD, IA-2, IAA long-term quality control for Diabetes type 1 autoantibodies

- use for internal quality control in medical laboratories
- 100% human serum
- 100% identical to patient sample
- available in low (L1) and high (L2) elevated range
- kit- and vendor independent
- shelf life up to 2 years
- production in line with highest ethical standards
- CE-marked
- volume: 1 ml per sample
- consignment of samples for testing possible



Price per sample 18.00 €  
(L1 (1 ml) or L2 (1 ml))



#### contact

**+49 3302-883 884** Mrs. Marleen Seiler  
**m.seiler@inventdiagnostica.de**  
**+49 3302-883 818** Mr. Ilja Strelkow  
**i.strelkow@inventdiagnostica.de**

Please give us noncommittal information

Please send us samples for testing

Our dates: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



certificate number:  
265396 QM

30.07.07